




LAB SOLUTION PROVIDER



TOPTION INSTRUMENT CO.,LTD



TOPTION GLASS




Workshop conner



RAW MATERIAL WAREHOUSE




GLASS COMPONENTS
DISPLAY ROOM





75°Bend type

Model	Grinding mouth
A151420	14/20
A151922	19/22
A152440	24/40



A type connector

Down - up, Large moth change to small

Model	Up grinding mouth	Down grinding mouth
A201419	14/20	19/22
A201424	14/20	24/40
A201924	19/22	24/40
A201929	19/22	29/42
A202429	24/40	29/42



B type connector

Down - up, small moth change to large

Model	Up grinding mouth	Down grinding mouth
A251914	19/22	14/20
A252414	24/40	14/20
A252419	24/40	19/22
A252914	29/42	14/20
A252924	29/42	24/40

Space-saving suction connector

Model	Grinding mouth	small joint outer diameter
A291420	14/20	8
A291922	19/22	8
A292440	24/40	10



Space-saving suction connector

#35/20 ball-milling connector with 10mm small joint at side

Model	Grinding mouth
A293520	35/20



Straight type suction connector

Straight grinding mouth connector with suction or gas guiding small joint at top

Model	Grinding mouth	small joint outer diameter
A301420	14/20	8
A301922	19/22	8
A302440	24/40	10





Straight type suction connector

Model	Grinding mouth	small joint outer diameter(mm)
A322440	24/40	10



Suction connector, double small joints

Top 180 ° direction has two small nozzle with opposite direction, used for nitrogen protection, could connect a bubbler.

Model	Grinding mouth	small joint outer diameter(mm)
A342440	24/40	10



Claisen connector

Model	Grinding mouth
A371420	14/20
A371922	19/22
A372440	24/40



90° Suction connector

PTFE valve connector with 90° small nozzle, outer diameter 8mm, inner hole diameter 2mm.

Model	Grinding mouth
A381420	14/20
A381922	19/22
A382440	24/40
A382942	29/42
Accessories	
S181125	2mm PTFE valve is optional



Lab Solution Provider



TOPTION INSTRUMENT CO.,LTD

Add:HeCheng, TaiBai Road, YanTa District, XI 'AN CHINA

Tel:0086-29-88763980

Fax:0086-29-88990306-8006

Email:info@toptionlab.com

info@toption-china.com

Web:www.toptionlab.com (for lab solution)

www.toption-china.com (for instrument)