

M series ultrasonic cleaner

M: machinery timing, no heating function

Model	Working size	Capacity	Frequency	Power	Timer	Cleaning basket	Lid	drainage
	L*W*H(mm)	L	KHz	W	min			
TP2-80M	150*140*100	2	40	80	1-30	√	√	N/A
TP3-120M	240*140*100	3	40	120	1-30	√	√	N/A
TP5-120M	240*140*150	5	40	120	1-30	√	√	√
TP4-180M	300*155*100	4	40	180	1-30	√	√	√
TP6-180M	300*155*150	6	40	180	1-30	√	√	√
TP10-250M	300*240*150	10	40	250	1-30	√	√	√
TP10-300M	300*240*150	10	40	300	1-30	√	√	√
TP14-400M	300*240*200	14.4	40	400	1-30	√	√	√
TP22-600M	500*300*150	22.5	40	600	1-30	√	√	√
TP30-800M	500*300*200	30	40	800	1-30	√	√	√

Ultrasonic cell crusher

Ultrasonic material dispenser

Ultrasonic Cleaner

