

## Household Mini Ultrasonic Cleaner(Plastic shell)

### Plastic shell M series ultrasonic cleaner

M: machinery timing, no heating function

Model	Working size	Capacity	Frequency	Power	Timer	Cleaning basket	lid	drainage
	L*W*H(mm)	L	KHz	W	min			
TCP-2M	200*138*75	2	40	60	1-99	√	√	N/A
TCP-3M	240*140*100	3	40	120	1-99	√	√	N/A
TCP-5M	240*140*150	5	40	120	1-99	√	√	√
TCP-10M	300*240*150	10	40	250	1-99	√	√	√



TCP-2M/A



TCP-3M/A



TCP-5M/A



TCP-10M/A

Ultrasonic cell crusher

Ultrasonic material dispenser

Ultrasonic Cleaner

### Plastic shell A series ultrasonic cleaner

A: single frequency, with heating function, LED display

Model	Working size	Capacity	Frequency	Power	Heating power	Temp range	Timer	Cleaning basket	Lid	drainage
	L * W * H (mm)	L	KHz	W	W	°C	min			
TCP-2A	200*138*75	2	40	60	100	RT-80	1-99	√	√	N/A
TCP-3A	240*140*100	3	40	120	100	RT-80	1-99	√	√	N/A
TCP-5A	240*140*150	5	40	120	120	RT-80	1-99	√	√	√
TCP-10A	300*240*150	10	40	250	300	RT-80	1-99	√	√	√