

Ultrasonic Cleaner industrial type

Product Description:

TP brand ultrasonic cleaner than with the cleaning of its own, but also with extraction, emulsification, dissolution speed, grinding, dispersion and other functions. Widely used: mechanical, electronics, plastics, instrumentation, environmental protection, medicine, packaging, military, aerospace, shipbuilding, automobile manufacturing and repair industry cleaning; experimental material straw, smoking Tsui and utensils clean before chromatography degassing treatment, medical equipment, medical materials and utensils clean; jewelry, jewelry, watches, precious metals, precious stones, coins, glasses cleaning



Ultrasonic cell crusher

Ultrasonic material dispenser

Ultrasonic Cleaner

Technical parameters:

Type	Ultrasonic Frequencies(kHz)	Ultrasonic Power (W)	Instrument Size (mm)	Working Size(mm)	Capacity (L)	Drainage	Lid	Timer	Heater
TP-1000T	28kHz	1000	620*520*550	500*400*300	60	√	optional	√	√
TP-1200T	28kHz	1200	620*520*550	500*400*300	60	√	optional	√	√
TP-1500T	28kHz	1500	770*600*550	650*480*300	90	√	optional	√	√
TP-1800T	28kHz	1800	820*600*550	700*480*300	100	√	optional	√	√
TP-2000T	28kHz	2000	850*580*560	750*480*300	110	√	optional	√	√
TP-2400T	28kHz	2400	900*700*560	800*600*300	145	√	optional	√	√

Optional frequency: 25kHz and 40kHz.