

C series ultrasonic cleaner
C: double frequency, with heating function, LCD display, ultrasonic power is adjustable

Model	Working size	Capacity	Frequency	Power	Power adjustable	Heating power	Temp range	Timer	Cleaning basket	Lid	drainage
	L * W * H (mm)	L	KHz	W	%	W	°C	min			
TP3-120C	240*140*100	3	25/45	120	1-100	150	RT-80	1-99	√	√	N/A
TP6-180C	300*155*150	6	25/45	180	1-100	180	RT-80	1-99	√	√	√
TP10-250C	300*240*150	10	25/45	250	1-100	300	RT-80	1-99	√	√	√
TP10-300C	300*240*150	10	25/45	300	1-100	300	RT-80	1-99	√	√	√
TP14-400C	300*240*200	14.4	25/45	400	1-100	300	RT-80	1-99	√	√	√
TP22-600C	500*300*150	22.5	25/45	600	1-100	500	RT-80	1-99	√	√	√
TP30-800C	500*300*200	30	25/45	800	1-100	700	RT-80	1-99	√	√	√

Ultrasonic cell crusher

Ultrasonic material dispenser

Ultrasonic Cleaner

