

## UP-4S-650W / UP-4S-1000W 4-channel ultrasonic crusher



### technical parameters:

Model	UP-4S-650W	UP-4S-1000W
Frequency(KHz)	20-25KHz Automatic tracking	
Power(W)	10-650W adjustable	10-1000W adjustable
Horn (mm)	Φ3*4 pcs	Φ6*4 pcs
Crush capacity at a time (ml)	3-15ml*4 kinds sample	10-150ml*4 kinds sample
Temperature alarm range	0-100 °C	
Instrument protection function	Over-temperature	
Foot switch	Optional	
Display	LCD display	