

A series ultrasonic cleaner
A: single frequency, with heating function, LED display

Model	Working size	Capacity	Frequency	Power	Heating power	Temp. range	Timer	Cleaning basket	Lid	drainage
	L * W * H (mm)	L	KHz	W	W	°C	min			
TP3-120A	240*140*100	3	40	120	150	RT-80	1-99	√	√	N/A
TP5-120A	240*140*150	5	40	120	150	RT-80	1-99	√	√	N/A
TP4-180A	300*155*100	4	40	180	180	RT-80	1-99	√	√	√
TP6-180A	300*155*150	6	40	180	180	RT-80	1-99	√	√	√
TP10-250A	300*240*150	10	40	250	300	RT-80	1-99	√	√	√
TP10-300A	300*240*150	10	40	300	300	RT-80	1-99	√	√	√
TP14-400A	300*240*200	14.4	40	400	300	RT-80	1-99	√	√	√
TP22-600A	500*300*150	22.5	40	600	500	RT-80	1-99	√	√	√
TP30-800A	500*300*200	30	40	800	700	RT-80	1-99	√	√	√

Ultrasonic cell crusher

Ultrasonic material dispenser

Ultrasonic Cleaner

